




Veterina Infusion Pump

New Design

- Exquisite and mini, po able
- 3.5"Large colo ul touch-screen and key-board, dual control methods.




Simple and Convenient

- Auto-lock Screen**
the lock screen duration can be adjusted in multiple levels; pre-vention for operation errors.
- Night Mode**
Design for night time safe operation.
- Histo Record Cassette**
Up to 50000pcs records.
Trackable treatments are ensured.
- System self-test**
System self-test is conducted eve time turning on, safeguard for eve -time use.

Accurate and Safe



- Large capacity lithium batte
- IP44 High-level dust-proof and waterproof



Con guration

Applicable Infusion Sets	1ml/15d, 1ml/19d, 1ml/20d, 1ml/60d
Flow Rate	ml/h mode: 0.1~1200ml/h (in 0.01ml/h increments) , drop/min mode: 1~(400ml/h*d/60) d/min (in 1d/min increments)
BOLUS	Bolus rate: 0.1~1200ml/h (in 0.01ml/h increments) , VTBI(Volume to be infused): 0.1~100ml (in 0.01ml increments)
FF (Fast Forward)	FF rate: 0.1~1200ml/h (in 0.01ml/h increments) , VTBI(Volume to be infused): 0~9999.99ml (in 0.01ml increments)
VTBI	0.1~9999.99ml (in 0.01ml increments)
Cumulative Infusion Volume	0~9999.99ml (in 0.01ml increments)
Rate (ml/h) Accuracy	±5% (after calibration);
Occlusion Pressure	5 occlusion pressure levels (20、50、80、110、130) kPa
Dimensions	131.5mm (L) ×101mm (W) ×163.8mm (H)