

Product name: B18 all stainless steel three-layer treatment cart

Size: 750*450*910mm



Parameters:

1. The frame is made of $\Phi 25 \times 1.1$ mm high-quality 304 stainless steel pipe, and the panel is made of 1.0mm high-quality 304 stainless steel brushed plate.
2. The upper, middle and lower layers can be used to place items for easy use, and the stainless steel solid core is fenced on three sides.
3. The top layer is equipped with a drawer, which uses a silent slide rail, which is easy and flexible to pull out and can store disposable infusion sets, cotton swabs, gauze, etc. Drawer top $\delta=1.0$ mm stainless steel brushed plate.

4. Four Taiwan Shi Ke Rui $\phi 100$ ultra-quiet casters, two of which have brakes for easy pushing.
5. Beautiful appearance, flat, straight, four parallel corners, no sharp edges, burrs and other obvious defects on the surface, all welded parts are polished flat and smooth, and polished evenly.
6. Suitable for placing instruments and items in the operating room.