



CM-11

Name: Drug Delivery Trolley

Model: CM-11

Specification: 750 * 475* 950mm

1. JEMP assembly line products are mainly composed of aluminum, steel, ABS engineering plastic structure; aluminum alloy triangular column four-column load-bearing;
2. The upper part of the car body: the large-scale mould of flying swallow type mesa including ABS guardrail gentle slope design is easy to clean, goods are not easy to slide down, the height of guardrail is 70mm, the use size in the mesa is 595 * 467mm, the use size in the mesa is 595 * 467mm, the use size in the mesa is 595 * 467mm, on the table with transparent soft glass, handrail guardrail dual-use integration design, after the implementation of arbitrary, double-layer design of the table on the right and left handrails, hollow handprint handle more beautiful and comfortable,
3. Front of the car body: center-locked rotary, five-layer medicine tray, three-fold light track, with a divider with plug-in label, a total of 25 four-box, each four-box can be divided into morning, middle, night, night. A total of 25 people, medicine box size: 216x100x84mm medicine box separate size: 100x48x84mm, a large draw 240mm, the inner space: 508 * 358 * 220mm, drawer 3 * 3 separators, can be freely separated, for all kinds of drugs and equipment storage;

4. To the left: side-hung lockers;
5. Right Side: flip color separation trash bin 175 * 175 * 280MM, double color for waste separation (yellow medical waste, green domestic waste)
6. On the back of the car body, there are 254-bit pill boxes with central control lock, 5-layer pill tray, 3-fold light sound track, and insert label with inner divider. The size of the pill box is 216 x 100 x 84 mm.
7. Body Bottom: four 100mm diameter luxury universal insert type Light Sound Wheel, two of which with brake function.