



CB-07

Based on in-depth research in the medical industry and considering the current complex medical environment, the medical R&D team has developed a highly practical and high-quality trolley: it can effectively help doctors conveniently and quickly transfer medical equipment, reduce unnecessary manpower consumption, and break down heavy and complex tasks.

Product Features

1. Mold-cast anti-slip edging, stable and sturdy, preventing machine slippage and falls.
2. High-quality aluminum alloy tabletop, easy to disinfect and clean.
3. Equipped with a large storage railing for convenient storage of

various items.

4. Die-cast aluminum alloy base, aesthetically pleasing and highly durable.

5. Rubber-coated casters, alkali-resistant, anti-slip, wear-resistant, smooth movement, and won't damage the floor.

6. Modular packaging, easy and convenient assembly.

Technical Specifications:

Load Capacity: 40kg

Trolley Height: 820mm

Trolley Weight: 10kg

Trolley Color: White

Packaging Dimensions: 190*860*670mm

Wheels: 7.5cm