



Specifications: 480*520*900/938/1138mm±10

Product Name: Laptop Mobile Nursing Cart/Ward Round Cart

1. Material: The car body is made of cold-rolled steel with powder coating. The upper table and drawers are made of ABS engineering plastic in one piece, which is wear-resistant and easy to clean. The aluminum alloy bottom is load-bearing, and it has a foot pedal for height adjustment.
2. Columns: Aluminum alloy outer column specifications 120*80mm, inner column 90*55mm, height adjustment range 210mm
3. Front: One-piece molded tabletop with double-layer armrests, ABS recessed guardrail 60mm. The first layer can hold a keyboard (460*200*50mm), left and right mouse worktables (200*190mm),

and the second layer is a custom battery box drawer.

4. Lower tabletop: 336*218*6mm

5. Left side: Left storage box with dividers, dimensions:
320*138*100/145mm.

6. Right side: Right wire basket, dimensions: 340*180*120/155mm,
comes with two 2-liter air purifiers.

7. Back: No accessories.

8. Bottom: Four-legged aluminum alloy mold, one-piece design,
ABS upper part, double-sided deluxe diameter silent wheels, two of
which have brakes, sturdy and durable; can be pushed and turned
freely on flat surfaces.

9. Packaging dimensions: 560*530*960mm

10. Gross weight: 34KG